Discovery Labware, Inc. 6 Henshaw Street Woburn, MA 01801 Tel: 1.978.442.2200 (U.S.) CLSTechServ@Corning.com

CORNING

www.corning.com/lifesciences

Insect Cell Control SUPERSOMES[™]

Catalog Number.......456200 Storage Conditions..STORE AT -80°C Date Released2014 May Expiration Date......2017 April

Insect cells (BTI-TN-5B1-4) infected with wild type baculovirus (*Autographa californica*) were used to prepare these microsomes.

ADVICE

- Thaw rapidly in a 37°C water bath. Keep on ice until use.
- Aliquot to minimize freeze-thawing cycles.
- When used as a control for native activities, use protein concentrations comparable to those for active microsome preparations.

Discovery Labware, Inc. 6 Henshaw Street Woburn, MA 01801 Tel: 1.978.442.2200 (U.S.) CLSTechServ@Corning.com

CORNING

www.corning.com/lifesciences

INSECT CELL MICROSOMES SAFETY INFORMATION

HAZARD WARNING:

The product was produced using baculovirus (*Autographa californica*) infected insect cells (BTI-TN-5B1-4). This virus is not known to be pathogenic to humans or other mammals.

SAFETY RECOMMENDATIONS:

When using this product, follow good laboratory safety procedures:

Do not eat, drink or smoke.

Avoid contact with skin or eyes.

Do not inhale aerosols.

Do not pipette by mouth.

Wear suitable protective clothing, gloves and eye protection.

Steam sterilize product or treat product with a 1% solution of sodium hypochlorite prior to disposal.