

## Corning® rLaminin-521 (Human)

PRODUCT: Corning rLaminin-521 (Human)  
CATALOG: 354221  
LOT NUMBER: 80154  
SOURCE: Recombinant protein expressed in a human cell line  
QUANTITY: 100 µg  
CONCENTRATION: 100 µg/ml  
FORMULATION: Frozen solution in Dulbecco's Phosphate-Buffered Saline (pH 7.2) with 10% glycerol, and 0.02% sodium azide.  
QUALITY CONTROL: rLaminin-521 (Human) is filtered through 0.2 micron sterile filter.

- **Authenticity:** Verified with Western Blot using specific antibodies against each subunit.
- **Purity:** ≥ 95% as determined by SDS-PAGE
- **Endotoxin Level:** ≤ 50 EU/ml
- **Mycoplasma:** Negative


STORAGE AND HANDLING: Store at -20°C. For short term use, rLaminin-521 can be stored at 2-8°C for 3-4 weeks. Make aliquots to avoid repeated freezing and thawing. Protein solution should be gently mixed before use.

EXPIRATION DATE: June 11, 2021, when stored at -20°C.

SAFETY RECOMMENDATION: Handle in accordance with good industrial hygiene and laboratory safety practices

REFERENCES:

1. Domogatskaya A., et al. Ann. Rev. Cell Dev. Biol. 28, 523 (2012)
2. Aumailley M., et al. Matrix Biol. 24, 326 (2005)
3. Rodin S., et al. Nature Comm. 5, 3195 (2014)

QUALITY ASSURANCE:  DATE: September 17, 2018

Discovery Labware, Inc., Two Oak Park, Bedford, MA 01730, Tel: 1.978.442.2200 (U.S.)  
CLSTechServ@Corning.com www.corning.com/lifesciences

CORNING

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

For a listing of trademarks, visit [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks)

© 2013 Corning Incorporated