

**Corning Incorporated
Life Sciences**

**Registered
ISO 9001:2008**

Certificate of Compliance

Catalog Number: 431231
Lot Number: 26212061
Product Description: Corning® syringe filter, Teflon™ (PTFE) membrane, 0.45 µM pore size, 25 mm diameter, bulk pack

Component Materials:

Housing - Polypropylene, meets *USP, Class VI* requirements for plastic containers and closures.
Membrane - Teflon™ (PTFE), 0.45 µM pore size, meets *USP, Class VI* requirements for plastic containers and closures.
Teflon is a registered trademark of E.I. DuPont DeNemours, Wilmington, Delaware

HPLC Certified:

The product has been tested to ensure there are no extraneous peaks in the UV range.

This manufacturing lot has been sampled and tested in accordance with Standard Operating Procedures and has been released by Quality Assurance for the following characteristics:

Test/Procedure:	Results:
Membrane Bubble Point Integrity	Pass
Housing Burst Pressure Test	Pass
Pressure Hold Test	Pass
Visual Attributes	Pass
Packaging	Pass

Expiration Date:

This product has an expiration date of 3 years from the date of manufacture for the attributes contained within this certificate.

Revision Date:

9-11-12



Wendy Loitsch
Quality Manager

Rev No: 2