

**Corning Incorporated  
Life Sciences**

**Registered  
ISO 9001:2008**

# Certificate of Compliance

**Catalog Number:** 431231  
**Lot Number:** 31912016  
**Product Description:** Corning® syringe filter, Teflon™ (PTFE) membrane, 0.45 µM pore size, 25 mm diameter, bulk pack

**Component Materials:**

Housing - Polypropylene, meets *USP, Class VI* requirements for plastic containers and closures.  
Membrane - Teflon™ (PTFE), 0.45 µM pore size, meets *USP, Class VI* requirements for plastic containers and closures.  
Teflon is a registered trademark of E.I. DuPont DeNemours, Wilmington, Delaware

**HPLC Certified:**

The product has been tested to ensure there are no extraneous peaks in the UV range.

This manufacturing lot has been sampled and tested in accordance with Standard Operating Procedures and has been released by Quality Assurance for the following characteristics:

<b>Test/Procedure:</b>	<b>Results:</b>
Membrane Bubble Point Integrity	Pass
Housing Burst Pressure Test	Pass
Pressure Hold Test	Pass
Visual Attributes	Pass
Packaging	Pass

**Expiration Date:**

This product has an expiration date of 3 years from the date of manufacture for the attributes contained within this certificate.

**Revision Date:**

9-11-12



Wendy Loitsch  
Quality Manager

Rev No: 2