

## Insect Cell Control SUPERSOMES™

Catalog Number.....456200  
Lot Number.....1253002

Storage Conditions..STORE AT -80°C  
Date Released .....2021 September  
Expiration Date.....2031 September

Package Contents.....1.5 mL  
Protein Content.....5.0 mg/mL in 100mM potassium phosphate (pH 7.4)  
Cytochrome c Reductase Activity.....30 nmole/(min x mg protein)  
Cytochrome P450 Content.....Not Detectable

Insect cells (BTI-TN-5B1-4) infected with wild type baculovirus (*Autographa californica*) were used to prepare these microsomes.

### ADVICE

- Thaw rapidly in a 37°C water bath. Keep on ice until use.
- Aliquot to minimize freeze-thawing cycles.
- When used as a control for native activities, use protein concentrations comparable to those for active microsome preparations.

***For research use only. Not for use in diagnostic or therapeutic procedures.***

For a listing of trademarks, visit [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks)

© 2013 Corning Incorporated


## INSECT CELL MICROSOMES

### HAZARD WARNING:

The product was produced using baculovirus (*Autographa californica*) infected insect cells (BTI-TN-5B1-4). This virus is not known to be pathogenic to humans or other mammals.

### SAFETY INFORMATION:

Safety assessment indicates this product is not hazardous, therefore no SDS (Safety Data Sheet) is provided. Use standard laboratory practices for the handling and disposal of Biosafety Level 1 materials.

  
\_\_\_\_\_  
Quality Assurance

13 Sept, 2021  
\_\_\_\_\_  
Date

***For research use only. Not for use in diagnostic or therapeutic procedures.***

For a listing of trademarks, visit [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks)

© 2013 Corning Incorporated