

Certificate of Compliance

Description: Filling Bell Assembly
2 ft. 1/4" I.D. Clear C-Flex Tubing w/ Female MPC
Connector w/ 1/4" Barb w/ Pinch Clamp and Filling Bell

Cat #: 91-700-12
Lot #: 44894-0000
Irradiation Date: 11/10
Expiration Date: 10/13

MSDS: Available upon request

Shipping: Ambient

Storage: Store at ambient temperature.

This is to certify that this product has been manufactured in a manner compliant with Current Good Manufacturing Practices (cGMP's), FDA Regulations 21 CFR part 820 and ISO 13485 Quality Standards. This product has been inspected and tested in accordance with the applicable material specifications and drawing requirements.

TEST PARAMETERS	SPECIFICATION	RESULT
Physical Inspection	Pass	Pass
Gamma Irradiation	Dosage: 27.5 kGy, Minimum 45.0 kGy, Maximum	Recorded Dosage: 28.8 kGy 40.4 kGy

NOTE: This product is intended for *Further Manufacturing/Processing* Use. Utilization of this product apart from the labeled intended use may be a violation of Federal Law. Mediatech, Inc. will not be liable under any warranty expressed or implied, in any amount greater than the price of the product sold.

Approved: *Seli Zurek* 05 NOV 2010
Quality Systems

Verified By/Date: *J. Levesy / Rauche* 05 NOV 2010

1.800.CELLGRO : TOLL FREE
703.471.5955 : LOCAL AREA
703.471.1086 : FACSIMILE

9345 DISCOVERY BOULEVARD MANASSAS, VIRGINIA