

Corning Incorporated - Life Sciences  
271 County Route 64  
Big Flats NY 14814 USA  
[www.corning.com/lifesciences](http://www.corning.com/lifesciences)  
Refer to website for regional contact information.

Page: 1 / 1

|                        |   |                         |              |
|------------------------|---|-------------------------|--------------|
| <b>Product Name</b>    | : AxyPrep Mag Plasmid Kit- Medium - 384 Preps |                         |              |
| <b>Catalog Number</b>  | : MAG-P-M                                     | <b>Manufacture Date</b> | : 2018-02-14 |
| <b>Lot ID</b>          | : 04518001                                    |                         |              |
| <b>Expiration Date</b> | : 2019-02-14                                  |                         |              |

**Quality Management System** - Complies with the current version of the ISO 9001 Standard.

**Quality Control Testing** - Representative production samples are collected and inspected in accordance with current applicable product specifications. Inspection records are reviewed and approved by qualified personnel for product release.

- This product met Corning Incorporated - Life Sciences' high standards of quality at the time of batch/lot release.



Elizabeth M (Betty) Hansen  
Quality Manager