

Corning Incorporated - Life Sciences  
271 County Route 64  
Big Flats NY 14814 USA  
[www.corning.com/lifesciences](http://www.corning.com/lifesciences)  
Refer to website for regional contact information.

Page: 1 / 1

<b>Product Name</b>	: AxyPrep Mag PCR Normalizer Kit- 50 mL - 384 Preps		
<b>Catalog Number</b>	: MAG-PCR-NM-50	<b>Manufacture Date</b>	: 2018-03-22
<b>Lot ID</b>	: 08118001		
<b>Expiration Date</b>	: 2019-03-22		

**Quality Management System** - Complies with the current version of the ISO 9001 Standard.

**Quality Control Testing** - Representative production samples are collected and inspected in accordance with current applicable product specifications. Inspection records are reviewed and approved by qualified personnel for product release.

- This product met Corning Incorporated - Life Sciences' high standards of quality at the time of batch/lot release.



Elizabeth M (Betty) Hansen  
Quality Manager