

836 North Street - Life Sciences
Tewksbury MA 01876 USA
www.corning.com/lifesciences
Refer to website for regional contact information.

Page: 1 / 1

| | | | |
|------------------------|---|-------------------------|--------------|
| Product Name | : AxyPrep Mag PCR Normalizer Kit- 50 mL - 384 Preps | | |
| Catalog Number | : MAG-PCR-NM-50 | Manufacture Date | : 2021-07-01 |
| Lot ID | : 18221002 | | |
| Expiration Date | : 2022-07-01 | | |

Quality Management System - Complies with the current version of the ISO 9001 Standard.

Quality Control Testing - Representative production samples are collected and inspected in accordance with current applicable product specifications. Inspection records are reviewed and approved by qualified personnel for product release.

- This product met Corning Incorporated - Life Sciences' high standards of quality at the time of batch/lot release.



ANDREW J (ANDY) FRANKS
CMO Quality Manager