

**Product description**

*Container 125mL PP WhiteScrewCap 38 0/CS*

**Catalogue number**

*TP52C-013*

**Lot number**

*761541*

**Manufacturing date**

*2020/09/10*

**Sterility Assurance Level**

Non sterile

**Expiration statement**

*N/A*

**Performance testing (AQL)**

*Inspection based on a statistically significant number of samples towards to internal or agreed specifications that determine Assurance Quality Level.*



Cédric MAES  
Plant Quality Manager, Borre  
maesc@corning.com